

PART OF THE **DURAN**[®] BOTTLE SYSTEM


ISOLATED

DURAN[®] DOUBLE WALLED BOTTLE



DWK LIFE SCIENCES
DURAN WHEATON KIMBLE

DURAN® GLS 80® DOUBLE WALLED BOTTLE WIDE MOUTH



The DURAN® double walled, wide mouth bottles GLS 80® incorporate an integral jacket that isolates the contents from the external environment. Heated or cooled liquids can be circulated through the jacket to control the temperature within the screw topped DURAN® bottle. The DURAN® double walled bottles offer a sealable and more flexible alternative to open topped jacketed beakers.

Cat. No.	Capacity (mL)	Thread	d (OD) (mm)	h (mm)	Pack Unit
with screw cap and pouring ring from PP (blue)					
24 256 44 51	500	GLS 80®	110	175	1
24 256 54 56	1 000	GLS 80®	110	275	1

